

5MP CMOS Microscope Camera

- Video, Image acquisition, labelling, contrast lighting & color adjusting, calibration, size changing, image rotation, sharpness control, etc
- Aluminium alloy body
- Frame Rate / Sensor Resolution
 - 5fps (2592 × 1944 Full Frame)
 - 21fps (1280 × 960)
 - 25fps (1024 × 768)
 - 38fps (640 × 480)
- Filter : RGB
- Sensor : Micron MT9P031 Size: 1/2.5", Pixel size : 2.2×2.2μm
- CMount, Eyepiece Adapters Included
- Computer Interface : USB2.0/480Mb/s
- Automatic/Manual Exposure Control & White Balance
- Includes 1.8 Meter USB Cable
- Measurement & Analysis Capabilities, 0.01mm Micrometer Included
- MAC OSX or PC - Win 2000/XP/Vista/Win 7/Win 8
- Weight: 220g

